


**THE MEDSTAR™ 200  
MONITORING SYSTEM**



**Simple monitoring system  
for blood pressure and weight**

**CYBERNET MEDICAL®**

A Revolution In Home Monitoring For Chronic Disease Patients



Digital Blood Pressure Monitor

Precision Health Scale

## MEDSTAR SPECIFICATIONS:


Size: 4.94" x 2.75" x 1.28"  
 Weight: 5.7 Ounces with Batteries  
 Battery: 2 AA Alkalines  
 Battery Life: 6 Months

## LCD DISPLAY:

Text: 2 Lines x 16 Characters Alphanumeric

## PHONE COMMUNICATIONS TO MEDSTAR SERVER

## OPTIONAL WEB-BASED DATA MANAGEMENT SOFTWARE



- Simple patient operation
- Electronic serial number for patient identification
- Durable, water resistant, and small in physical size
- Independent battery power
- Stores data from device reading until successfully transferred
- Records time and date of all measurements

The MedStar 200 is currently optimized for use with the following measurement devices. Other devices may be supported. Please see [www.cybernetmedical.com](http://www.cybernetmedical.com) for details.\*

	Manufacturer	Model	Outputs
Non-Invasive Blood Pressure (NIBP)	A&D Medical	UA-767PC	Systolic Pressure (mmHg) Diastolic Pressure (mmHg) Pulse Rate (per minute)
Weight Scale	A&D Medical	UC-321 PL	Weight (up 450 lbs. or 205 kg.)

\* Ask us if you would like support for alternative devices.