


**THE MEDSTAR™ INR
& TEMP SERIES**



**Simple monitoring system
for INR and thermometer devices**

MEDStar™ INR and TEMP Series



MEDSTAR SPECIFICATIONS:

Size: 4.94" x 2.75" x 1.28"

Weight: 5.7 Ounces with Batteries

Battery: 2 AA Alkalines


Battery Life: 6 Months

LCD DISPLAY:

Text: 2 Lines x 16 Characters Alphanumeric

PHONE COMMUNICATIONS TO MEDSTAR SERVER

OPTIONAL WEB-BASED DATA MANAGEMENT SOFTWARE



- Simple patient operation
- Electronic serial number for patient identification
- Durable, water resistant, and small in physical size
- Independent battery power
- Stores data from device reading until successfully transferred
- Records time and date of all measurements

The MedStar 220 is currently optimized for use with the following measurement devices. Other devices may be supported. Please see www.cybernetmedical.com for details.*

	Manufacturer	Model	Outputs
International Normalized Ratio (INR) Meter	HemoSense	INRatio	Prothrombin time and INR results
Sure Temp Plus	Welch Allyn	690	Temperature (F° and C°)

INR/Temp Model Configurations			
220	230	360	370
NIBP INR Scale Thermometer	INR Spirometer Thermometer PulseOx	NIBP Thermometer INR Glucose	NIBP INR Scale Glucose

*Ask us if you would like support for alternative devices.